

- Asepro 560 disposable sterile protective hooded coverall can be safely used in sterile production areas, the pharmaceutical industry, and laboratories.
- It provides reliable protection against infectious agents, contamination, and cytostatic drugs.
- Its breathable fabric and ergonomic design ensure comfort even during long periods of use.

ASP560 HCA [Ⓟ] [Ⓒ] [Ⓢ] ¹

Protective Hooded Coverall



Elasticated hood that fully encloses facial equipment.



Elasticated wrists.



Elasticated ankles.



Product Properties

Hazards	Infective agents, Contamination, Cytostatic drugs
Controlled Environments	ISO Class 5 GMP A/B, Sterile
Materials	Laminated SS (55 g/m ²)
Seams	Serged
Sizes	S, M, L, XL, 2XL, 3XL, 4XL ²
Colors	●
Packing	SD: Sterile, individually double packed
Quantity per Bag	1
Quantity per Box	50

Notes

¹ Ⓟ : Packing, Ⓒ : Color, Ⓢ : Size.
Learn more about the part numbering system: [PNS-CL-en.pdf](#)
² Check the size table: [SZT-CL-en.pdf](#)

Asepro

Superior comfort, reliable protection, and cleanroom compatibility for controlled environments.



Comfort

All-day comfort and protection against dust and liquid splashes with a flexible and breathable structure.



Chemguard

High-level industrial protection against chemical liquids and gases.



www.coverlab.com.tr

ANECO