

- Asepro 600 disposable sterile protective sleeves can be safely used in sterile production areas, the pharmaceutical industry, and laboratories.
- It provides reliable protection against infectious agents, contamination, and cytostatic drugs.
- Its breathable fabric and ergonomic design ensure comfort even during long periods of use.

## Product Properties

Hazards	Infective agents, Contamination, Cytostatic drugs
Controlled Environments	ISO Class 5 GMP A/B, Sterile
Materials	Laminated SS (55 g/m <sup>2</sup> )
Seams	Serged
Sizes	UNI <sup>2</sup>
Colors	●
Packing	SD: Sterile, individually double packed
Quantity per Bag	1
Quantity per Box	200

### Notes

<sup>1</sup> ⊕ : Packing, ⊙ : Color, ⊗ : Size.  
Learn more about the part numbering system: [PNS-CL-en.pdf](#)  
<sup>2</sup> Check the size table: [SZT-CL-en.pdf](#)

For informational purposes only.  
Not yet available for sale.



CATEGORY III

CE 2841



TYPE PB [6-B]



EN 1073-2



EN 14126



ISO 11137



Sterile



Elasticated wrists.

# Asepro 600

## ASP600 SLV ⊕ ⊙ ⊗<sup>1</sup> Protective Sleeves

## Asepro

Superior comfort, reliable protection, and cleanroom compatibility for controlled environments.



## Comfort

All-day comfort and protection against dust and liquid splashes with a flexible and breathable structure.



## Chemguard

High-level industrial protection against chemical liquids and gases.

[www.coverlab.com.tr](http://www.coverlab.com.tr)

ANECO